

# TAbic® V.H. Live Newcastle Vaccine



Safe & Proven protection in a convenient presentation

TAbic® family - TAbic® V.H., TAbic® IBVAR206 and TAbic® IB Var, TAbic® H-120 and TAbic® M.B.

www.pahc.com

Code # GV18006/INT/0418





## TAbic<sup>®</sup> V.H.

TAbic® V.H. is a live virus vaccine for Newcastle disease, developed from the cloned V.H. strain, freeze-dried and tableted in an effervescent form, sealed in eco-friendly sterile blister packaging.

## Efficacy Proven protection against virulent ND field viruses all over the world Highly immunogenic - produces

## Safety

- Mild post-vaccination reactions
- Suitable for Broilers, Layers, Breeders and Turkeys at any age


#### Convenience

- The tablets dissolve rapidly and completely, resulting in an even distribution of the virus in the diluent
- Applicable by aerosol, coarse spray and drinking water

## **Efficacy**

high titers


TAbic<sup>®</sup> V.H. vaccinated SPF chicks (at 1 & 17 days of age) by spray and by oral routes, had higher antibody titers and higher survivability post challenge than chicks vaccinated with an enterotropic vaccine (South African trial, 2004)



Challenge Strain: Local SA no. HOD 45/95 ICPI = 1.85 genotype VI Velogenic viscerotropic (IM route)

## Safety

The V.H. strain is a very mild virus, as seen in the ICPI test, yet provides good protection



Test results in PAHC files

## TAbic® V.H. Advantages

#### Mild, safe and immunogenic

Mild post-vaccination reaction, highly immunogenic vaccine. Reliable protection against vvNDV challenges can be applied at any age

#### Innovative

Eco-friendly, freeze-dried effervescent tablets, sealed in aluminum blister packs

#### Easy to use

The tablets dissolve and rapidly disperse the vaccine into a homogeneous clear solution

#### User friendly

TAbic® V.H. is user friendly, convenient and safe to handle. No need for working tools, syringes or other equipment

#### **Environmentally friendly**

TAbic® V.H. is environmentally friendly, empty blisters are easy to dispose of

#### Saving cost on shipment and storage

Due to reduced pack size, TAbic® V.H. requires minimal storage space, is lightweight and easily transported

#### Storage

TAbic® V.H. storage temperature should be 2°C to 8°C



